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Dear Sir or Madam, thank you for purchasing our product!

Before using the product, please read the following instructions to ensure proper use of the product.

Please keep this manual for future reference and follow its recommendations, as failure to follow them may pose a risk to life or health.

APPLICATION AND DESCRIPTION OF THE DEVICE

This massage mat is designed for massaging the back, shoulders, hips, and thighs. With 5 built-in massage programs , 9 vibration points , and a heating function , the device helps reduce muscle tension, improves blood circulation, and promotes relaxation. It can be used at home, in the office, or in the car.

This device is intended for home use only and is not intended for medical or commercial use, and is not intended to replace medical therapy. The device is intended for indoor use in dry rooms. It is not intended for medical, therapeutic, or commercial use.

For safety and CE certification reasons, this product must not be rebuilt or modified in any way. Using the product for purposes other than those described may result in damage to the product. Improper use may also result in hazards such as short circuits, fire, electric shock, etc.

WARNING

Danger of improper use!

Using the device in a manner other than its intended purpose and/or using the device in any other way may involve various risks.

The device should only be used for its intended purpose.

Please follow the procedures described in this instruction manual.

Note : It is recommended to consult a doctor before using the massager, especially if you have inflammation, joint injuries or skin problems.

START-UP/INSTALLATION INSTRUCTIONS

- Spread the mat on an armchair, sofa or other stable surface.
- Plug the AC adapter into an outlet or use a 12V car adapter.
- Plug the power supply cord into the socket on the mat.
- Turn on the device using the PWR button on the remote control.
- Select the desired massage program (MODE button) and zone (shoulders, back, hips, thighs)
- Set the intensity (INT button).
- If necessary, activate the heating function using the HEAT button .
- Set the operating time using the TIMER button. The mat will automatically turn off after the cycle is complete.

TIP

- Using the device in a manner other than its intended purpose and/or using the device in any other way may involve various risks.
- Use the device only for its intended purposes.
- The procedures described in this operating manual must be followed. Claims of any kind resulting from damage caused by improper use, improper repairs, unauthorized modifications, or the use of unauthorized spare parts are excluded. The risk lies solely with the user.
- If you experience any discomfort or skin irritation, discontinue use immediately.

TECHNICAL DATA

- Heating function: YES
- Sleep timer: YES (15,30,60 minutes)
- 5 massage programs
- 9 massage points
- 4 massage zones
- 3 levels of intensity
- Remote control included
- Power supply: mains AC 100V – 240V 50/60Hz + car

SAFETY INSTRUCTIONS

- This product may be used by children aged 8 years and over and by persons with reduced physical, sensory or mental capabilities or lack of experience or knowledge, provided they are supervised or have been instructed on the safe use of the product and understand the resulting risks.
- Never immerse the device in water.
- For cleaning, use a damp cloth or a mild detergent.
- Do not use a damaged device.

- The product is not intended for use by children.
- Children must not play with the appliance.
- Cleaning and user maintenance must not be performed by children without supervision.
- The product itself and its accessories are not toys and should be kept out of the reach of children to avoid danger.
- Make sure packaging materials are not left unattended. Children may play with them, which is dangerous.
- Protect the product from extreme temperatures, direct sunlight, strong vibrations, high humidity, moisture, flammable gases, vapors and solvents.
- Do not expose the product to mechanical stress.
- If safe operation is no longer possible, discontinue use and secure the product against further use. Safe operation is no longer possible if the product: - has been damaged, - is not functioning properly, - has been stored for an extended period under unfavorable conditions, or - has been subjected to excessive stress during transport.
- Do not use the product if any part is damaged. If the cord is damaged, do not attempt to repair it yourself.
- Do not use if the power cord or power adapter is damaged or if it is not properly plugged into the outlet.
- Do not disassemble the device yourself.
- Do not use on open wounds, bruises, swelling, or skin lesions. This may worsen the condition.
- Do not use in a vehicle while driving. The mat should only be used in a vehicle when it is stationary.
- Unplug the device after each use. Do not leave it plugged in while unattended.
- Do not fall asleep while the device is operating. This may result in excessive pressure and overheating of the body.
- Do not use the device immediately after a meal. It is recommended to wait at least 30 minutes.
- Do not use the device after consuming alcohol, sleeping pills, or sedatives. This may impair your response and increase the risk of injury.
- Discontinue the massage if you experience pain, dizziness, or discomfort. This device is intended for relaxation, not medical therapy.
- Before using the device, make sure that the mains voltage matches the required operating voltage of the device.
- The product should always be used as intended.

ELECTRIC SHOCK HAZARD

- Connect the device only to a properly installed and grounded power outlet. The mains voltage must match the data on the device/power supply rating plate.
- Ensure that the power cord does not become wet or damp during operation. Route the cord so that it does not become pinched or damaged.
- Keep the power cable away from hot surfaces.
- Do not perform any repairs to the device yourself.
- Any repairs may only be performed by service or qualified specialists.
- Before cleaning the device, unplug the power supply from the mains socket.
- Do not use the device in a damp or wet environment, such as a bathroom.
- Have a damaged plug or power cable repaired immediately by a qualified specialist or service center to avoid any hazards.
- Do not insert any needles or sharp objects into the device.

WARNING! RISK OF INJURY!

- This device should not be used to massage infants, small children, helpless people or animals.
- Never use this device if you experience dislocation or pain in the area being massaged.
- Stop the massage immediately if it becomes uncomfortable or causes pain.
- The device has a hot surface during operation. People who are not sensitive to heat should exercise caution when using the device.
- Do not use the device on swollen, burned, irritated, diseased or injured parts of the skin or body.
- If you have any doubts regarding use, please consult your doctor.
- Do not cover the device (pillow, blanket, etc.) during use.

WARNING – MATERIAL DAMAGE!

- The device is for indoor use only.
- The device should only be used with the supplied mains adapter.

TIPS AND INFORMATION ON THE MANAGEMENT OF USED PACKAGING

The packaging is made of environmentally friendly materials that can be disposed of at your local recycling center.

Used packaging material should be delivered to a designated waste disposal site designated by local authorities. For information on recycling options for used products, please contact your local municipality or city office.

DISPOSAL OF USED ELECTRICAL AND ELECTRONIC EQUIPMENT

For environmental reasons, used electrical and electronic products should not be disposed of with regular household waste, but should be disposed of properly. Information on collection points and their opening hours can be obtained from the relevant authorities.

This product complies with the requirements of the applicable European and national directives. The product meets European and national requirements for the safety of devices and products.

This product complies with the requirements of the relevant European and national RoHS directives.

Warning of electric shock! Danger to life!

We reserve the right to make changes to the text, design and technical data of the product without notice.

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